

KIT20P

The KIT20P is a *Cellbag*[®] Holder kit for the *Wave Bioreactor*[®] *BASE20/50P*. The KIT20P includes a white mounting plate, TOP-PLATE20, and a stainless steel tray, the HOLDER20. The HOLDER20 is used to support the flexible bioreactor chamber during transportation and operation. It is a non-disposable component. The *Cellbag* disposable bioreactor locks into the holder which in turn rests on the rocker unit (BASE20/50P). The holder conforms the inflated *Cellbag* to the optimal geometry for oxygen transfer and mixing. It also provides a containment vessel in the unlikely event of a leak.

Either two 2 liter *Cellbags*, two 10 liter *Cellbags*, or a single 20 liter *Cellbag* may be used on a KIT20P.

A protective PVC cover (LID20 - clear or LID20W - opaque) is available for benchtop operation.

Performance:	High quality polished stainless steel. Easy to use locking assembly for Cellbags Leak containment
Dimensions:	711mm x 575mm x 77mm (28" x 22.6" x 3")
Weight:	HOLDER20 5.4Kg (12 lb) KIT20P 9.3Kg (20.5 lb)
Power:	None