

## LID20W Cellbag<sup>®</sup> Cover (Opaque)

The LID20W is an opaque white polyvinyl chloride (PVC) cover for the KIT20EH, KIT20EHT, KIT20EHTD and KIT20P used with the *Wave Bioreactor<sup>®</sup> System20/50 series*. It is used to protect the flexible bioreactor chamber during benchtop operation and fits over the HOLDER20 and HOLDER20T. It is made in this opaque white version for light-sensitive applications.



<b>Performance:</b>	High quality opaque PVC construction Access hatch for sampling and additions May be cleaned with isopropyl alcohol or bleach
<b>Dimensions:</b>	725mm (L) x 580mm (D) x 210mm (H) (28.5" x 22.8" x 8.25")
<b>Weight:</b>	2.9 Kg (6.4 lbs)
<b>Power:</b>	None
<b>Options:</b>	None