

LID20 Cellbag[®] Cover

The LID20 is a clear polyvinyl chloride (PVC) cover for the KIT20EH, KIT20EHT, KIT20EHTD and KIT20P used with the *Wave Bioreactor[®] System20/50 series*. It is used to protect the flexible bioreactor chamber during benchtop operation and fits over the HOLDER20 or HOLDER20T.



Performance: High quality clear PVC construction
Access hatch for sampling and additions
May be cleaned with isopropyl alcohol or bleach

Dimensions: 725mm (L) x 580mm (D) x 210mm (H) (28.5" x 22.8" x 8.25")

Weight: 2.9 Kg (6.4 lbs)

Power: None

Options: None