

LID50W Cellbag[®] Cover (Opaque)

The LID50W is an opaque white polyvinyl chloride (PVC) cover for the KIT50EH, KIT50EHT, KIT50EHTD and KIT50P used with the *Wave Bioreactor[®] System20/50 series*. It is used to protect the flexible bioreactor chamber during benchtop operation and fits over the HOLDER50 and HOLDER50T. It is made in this opaque white version for light-sensitive applications.



Performance: High quality opaque PVC construction
Access hatch for sampling and additions
May be cleaned with isopropyl alcohol or bleach

Dimensions: 790mm (L) x 720mm (D) x 225mm (H) [31" x 28.3" x 8.8"]

Weight: 4.5 Kg (9.9 lbs)

Power: None

Options: None