

MIXHOLDER20 Bag Holder

MIXHOLDER20 is a fiberglass bag holder for Flexboys® and other disposable plastic bags for use on the Wave Mixer. It is used to support the flexible bags during transportation and operation. It is a non-disposable component. The bags are simply placed into the holder which in turn rests on to the rocker unit (MIXER20/50). The holder conforms the bag to the optimal geometry for mixing. It also provides a containment vessel in the unlikely event of a leak.



Performance:	High quality NSF certified molded white fiberglass. Leak containment.
Dimensions:	660mm x 460mm x 75mm
Weight:	2Kg (4.6lb)
Power:	None.
Options:	Fiberglass Lid (MIXLID20).