

## MIXHOLDER50 Bag Holder

A fiberglass bag holder for Flexboys® and other disposable plastic bags must be used on the Wave Mixer. It is used to support the flexible bags during transportation and operation. It is a non-disposable component. The bags are simply placed into the holder which in turn rests on to the rocker unit (MIXER20/50). The holder conforms the bag to the optimal geometry for mixing. It also provides a containment vessel in the unlikely event of a leak.



<b>Performance:</b>	High quality NSF certified molded white fiberglass. Leak containment.
<b>Dimensions:</b>	740mm x 635mm x 128mm
<b>Weight:</b>	4Kg (8.8lb)
<b>Power:</b>	None.
<b>Options:</b>	None