

## MIXKIT1000E

Stainless-steel bag holder for M\*Bag 1000 disposable plastic mixing bags. For use with Wave Mixer500/1000E base rocker unit. It is used to support the flexible bags during mixing and transportation. It is a non-disposable component. The bags are simply placed into the holder which in turn rests on to the rocker unit. The holder conforms the bag to the optimal geometry for mixing. It also provides a containment vessel in the unlikely event of a leak.

<b>Performance:</b>	High quality all stainless steel construction. Quick-release bag holddown clamps. Sealed construction for leak containment.
<b>Dimensions:</b> (w x d x h)	61 in x 91 in x 22 in (155 cm x 231 cm x 56 cm)
<b>Weight:</b>	750 lbs (340 Kg)
<b>Power:</b>	4 x 1250W @ 230 VAC for optional heaters
<b>Options:</b>	Heaterpads (Select MIXKIT1000EH for use with MIXER500/1000EH).