

STERILE TUBE FUSER – PRINTER

The STF-PRINTER can be connected to the printer port on the STERILE TUBE FUSER-IRc. It will print the weld settings, tube type/size selected, and status of the tube fusing operation. The data and time of each operation is also printed. If an error or failure occurs during the weld cycle it will be printed. The STF-PRINTER is an ideal accessory for validation and GMP documentation.



```

=====
UNIT ID: 0025      WELD # 00021

PROG ver: 2,1      DATA ver: 2,0
TUBE TYPE: CFLEX  5/8 in OD
----- SETTINGS -----
STER TEMP      400 C
STER TIME      10,0 s
CUT TEMP       220 C
WELD TEMP      240 C
COOL TEMP      50 C
CUT DLY        1,0 s
ALIGN DLY      0,0 s
WELD TIME      20,0 s
RELS TIME      5,0 s
XMOVE 1        0 counts
XMOVE 2        200 counts
XMOVE 3        1200 counts
XTIME 3        1,9 s
----- PROCESS -----
WELD START     11/01/2004 10:43
STERILIZED     11/01/2004 10:43
WELD COMPLETE  11/01/2004 10:45
=====

```

SAMPLE PRINTOUT

- Catalog Number:** STF-IRc.003214
- Performance:** Thermal 32 column printer
Serial RS232 interface
Auto paper load
Serial Printer cable included
- Dimensions:** 145mm W x 195mm D x 121mm H (5.7" x 7.7" x 4.8")
- Weight:** 1.2 Kg (2.7 lbs)
- Power:** 24 VDC, 110/220 VAC auto ranging adapter included.
- Options:** Reorder paper roll STF-IRc.003215 (6 rolls/pkg.)